

Product Lineup

Material: Fluoropolymer
Recommended Temperature: 200°C / 392° F

NFPX		20:1 Shrink Ratio / Standard / Standard Length: 1000 mm / Color: Transparent					MOQ/SPQ:100/100 [pcs]
Parts No.	Metric			Imperial			French
	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Suitable for
NFPX0045	0.60 mm	0.30 mm	0.20 +/- 0.05 mm	0.0236"	0.0118"	0.0079 +/- 0.0020"	1.3 - 1.5 Fr
NFPX0090	1.20 mm	0.60 mm	0.25 +/- 0.05 mm	0.0472"	0.0236"	0.0098 +/- 0.0020"	2.2 - 2.3 Fr
NFPX0150	2.00 mm	1.00 mm	0.275 +/- 0.05 mm	0.0787"	0.0394"	0.0108 +/- 0.0020"	3.4 - 5.7 Fr

NFPR		1.8:1 Shrink Ratio / Standard / Standard Length: 1700 mm / Color: Translucent					MOQ/SPQ:100/100 [pcs]
Parts No.	Metric			Imperial			French
	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Suitable for
NFPR006	0.72 mm	0.40 mm	0.25 +/- 0.05 mm	0.0283"	0.0157"	0.0098 +/- 0.0020"	1.6 - 1.8 Fr
NFPR007	0.90 mm	0.50 mm	0.25 +/- 0.05 mm	0.0354"	0.0197"	0.0098 +/- 0.0020"	1.9 - 2.4 Fr
NFPR008	1.08 mm	0.60 mm	0.25 +/- 0.05 mm	0.0425"	0.0236"	0.0098 +/- 0.0020"	2.2 - 2.9 Fr
NFPR010	1.26 mm	0.70 mm	0.25 +/- 0.05 mm	0.0496"	0.0276"	0.0098 +/- 0.0020"	2.5 - 3.5 Fr
NFPR011	1.44 mm	0.80 mm	0.25 +/- 0.05 mm	0.0567"	0.0315"	0.0098 +/- 0.0020"	2.8 - 4.0 Fr
NFPR013	1.62 mm	0.90 mm	0.25 +/- 0.05 mm	0.0638"	0.0354"	0.0098 +/- 0.0020"	3.1 - 4.6 Fr
NFPR014	1.80 mm	1.00 mm	0.25 +/- 0.05 mm	0.0709"	0.0394"	0.0098 +/- 0.0020"	3.4 - 5.1 Fr
NFPR018	2.34 mm	1.30 mm	0.25 +/- 0.05 mm	0.0921"	0.0512"	0.0098 +/- 0.0020"	4.3 - 6.7 Fr
NFPR022	2.88 mm	1.60 mm	0.28 +/- 0.05 mm	0.1134"	0.0630"	0.0110 +/- 0.0020"	5.2 - 8.4 Fr
NFPR029	3.78 mm	2.10 mm	0.30 +/- 0.05 mm	0.1488"	0.0827'	0.0118 +/- 0.0020"	6.7 - 11.1 Fr

NFPM		1.6:1 Shrink Ratio / Ultra Small / Standard Length: 1000 mm / Color: Translucent					MOQ/SPQ:100/100 [pcs]
Parts No.	Metric			Imperial			French
	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Suitable for
NFPM0030	0.38 mm	0.23 mm	0.20 +/- 0.05 mm	0.0150"	0.0091"	0.0079 +/- 0.0020"	0.9 - 1.0 Fr
NFPM0037	0.46 mm	0.28 mm	0.20 +/- 0.05 mm	0.0181"	0.0110"	0.0079 +/- 0.0020"	1.0 - 1.2 Fr
NFPM0042	0.51 mm	0.33 mm	0.20 +/- 0.05 mm	0.0201"	0.0130"	0.0079 +/- 0.0020"	1.2 - 1.3 Fr

NFPSC		1.6:1 Shrink Ratio / Standard / Length: 1700 mm / Color: Transparent					MOQ/SPQ:100/100 [pcs]
Parts No.	Metric			Imperial			French
	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Suitable for
NFPSC005	0.64 mm	0.40 mm	0.20 +/- 0.05 mm	0.0252"	0.0157"	0.0079 +/- 0.0020"	1.6 Fr
NFPSC007	0.80 mm	0.50 mm	0.20 +/- 0.05 mm	0.0315"	0.0197"	0.0079 +/- 0.0020"	1.9 - 2.1 Fr
NFPSC008	0.96 mm	0.60 mm	0.20 +/- 0.05 mm	0.0378"	0.0236"	0.0079 +/- 0.0020"	2.2 - 2.6 Fr
NFPSC009	1.12 mm	0.70 mm	0.20 +/- 0.05 mm	0.0441"	0.0276"	0.0079 +/- 0.0020"	2.5 - 3.0 Fr
NFPSC010	1.28 mm	0.80 mm	0.20 +/- 0.05 mm	0.0504"	0.0315"	0.0079 +/- 0.0020"	2.8 - 3.5 Fr
NFPSC012	1.44 mm	0.90 mm	0.20 +/- 0.05 mm	0.0567"	0.0354"	0.0079 +/- 0.0020"	3.1 - 4.0 Fr
NFPSC013	1.60 mm	1.00 mm	0.20 +/- 0.05 mm	0.0630"	0.0394"	0.0079 +/- 0.0020"	3.4 - 4.5 Fr
NFPSC017	2.08 mm	1.30 mm	0.20 +/- 0.05 mm	0.0819"	0.0512"	0.0079 +/- 0.0020"	4.3 - 6.0 Fr
NFPSC021	2.56 mm	1.60 mm	0.275 +/- 0.05 mm	0.1008"	0.0630"	0.0108 +/- 0.0020"	5.2 - 7.4 Fr
NFPSC028	3.36 mm	2.10 mm	0.275 +/- 0.05 mm	0.1323"	0.0827"	0.0108 +/- 0.0020"	6.7 - 9.8 Fr

NFPS		1.6:1 Shrink Ratio / Standard / Standard Length: 1700 mm / Color: Translucent					MOQ/SPQ:100/100 [pcs]
Parts No.	Metric			Imperial			French
	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Suitable for
NFPS005	0.64 mm	0.40 mm	0.20 +/- 0.05 mm	0.0252"	0.0157"	0.0079 +/- 0.0020"	1.6 Fr
NFPS007	0.80 mm	0.50 mm	0.20 +/- 0.05 mm	0.0315"	0.0197"	0.0079 +/- 0.0020"	1.9 - 2.1 Fr
NFPS008	0.96 mm	0.60 mm	0.20 +/- 0.05 mm	0.0378"	0.0236"	0.0079 +/- 0.0020"	2.2 - 2.6 Fr
NFPS009	1.12 mm	0.70 mm	0.20 +/- 0.05 mm	0.0441"	0.0276"	0.0079 +/- 0.0020"	2.5 - 3.0 Fr
NFPS010	1.28 mm	0.80 mm	0.20 +/- 0.05 mm	0.0504"	0.0315"	0.0079 +/- 0.0020"	2.8 - 3.5 Fr
NFPS012	1.44 mm	0.90 mm	0.20 +/- 0.05 mm	0.0567"	0.0354"	0.0079 +/- 0.0020"	3.1 - 4.0 Fr
NFPS013	1.60 mm	1.00 mm	0.20 +/- 0.05 mm	0.0630"	0.0394"	0.0079 +/- 0.0020"	3.4 - 4.5 Fr
NFPS017	2.08 mm	1.30 mm	0.20 +/- 0.05 mm	0.0819"	0.0512"	0.0079 +/- 0.0020"	4.3 - 6.0 Fr
NFPS021	2.56 mm	1.60 mm	0.275 +/- 0.05 mm	0.1008"	0.0630"	0.0108 +/- 0.0020"	5.2 - 7.4 Fr
NFPS028	3.36 mm	2.10 mm	0.275 +/- 0.05 mm	0.1323"	0.0827"	0.0108 +/- 0.0020"	6.7 - 9.8 Fr

NFPS		1.6~1.4:1 Shrink Ratio / Large / Standard Length: 1000 mm / Color: Translucent					MOQ/SPQ:50/50 [pcs]
Parts No.	Metric			Imperial			French
	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Suitable for
NFPS036	4.45 mm	2.79 mm	0.275 +/- 0.05 mm	0.1752"	0.1098"	0.0108 +/- 0.0020"	8.8 - 13.1 Fr
NFPS038	4.78 mm	2.92 mm	0.275 +/- 0.05 mm	0.1882"	0.1150"	0.0108 +/- 0.0020"	9.2 - 14.1 Fr
NFPS052	6.05 mm	4.32 mm	0.275 +/- 0.05 mm	0.2382"	0.1701"	0.0108 +/- 0.0020"	13.4 - 18.0 Fr

NFP		1.5:1 Shrink Ratio / Extra Large / Standard Length: 1500 mm / Color: Translucent					MOQ/SPQ:25/25 [pcs]
Parts No.	Metric			Imperial			French
	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Suitable for
NFP064	7.62 mm	5.08 mm	0.33 +/- 0.05 mm	0.3000"	0.2000"	0.0130 +/- 0.0020"	15.7 - 22.7 Fr
NFP074	8.89 mm	5.93 mm	0.33 +/- 0.05 mm	0.3500"	0.2335"	0.0130 +/- 0.0020"	18.3 - 26.6 Fr
NFP085	10.16 mm	6.77 mm	0.33 +/- 0.05 mm	0.4000"	0.2665"	0.0130 +/- 0.0020"	20.9 - 30.4 Fr
NFP095	11.43 mm	7.62 mm	0.33 +/- 0.05 mm	0.4500"	0.3000"	0.0130 +/- 0.0020"	23.4 - 34.3 Fr

NFP		1.4:1 Shrink Ratio / Standard / Standard Length: 1500 mm / Color: Translucent					MOQ/SPQ:100/100 [pcs]
Parts No.	Metric			Imperial			French
	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Expanded I.D. Min.	Full Recovered I.D. Max.	Full Recovered Wall Thickness	Suitable for
NFP006	0.71 mm	0.50 mm	0.20 +/- 0.05 mm	0.0280"	0.0197"	0.0079 +/- 0.0020"	1.8 - 1.9 Fr
NFP007	0.86 mm	0.60 mm	0.20 +/- 0.05 mm	0.0339"	0.0236"	0.0079 +/- 0.0020"	2.2 - 2.3 Fr
NFP008	1.00 mm	0.70 mm	0.20 +/- 0.05 mm	0.0394"	0.0276"	0.0079 +/- 0.0020"	2.5 - 2.7 Fr
NFP010	1.14 mm	0.80 mm	0.20 +/- 0.05 mm	0.0449"	0.0315"	0.0079 +/- 0.0020"	2.8 - 3.1 Fr
NFP011	1.29 mm	0.90 mm	0.20 +/- 0.05 mm	0.0508"	0.0354"	0.0079 +/- 0.0020"	3.1 - 3.6 Fr
NFP013	1.57 mm	1.10 mm	0.20 +/- 0.05 mm	0.0618"	0.0433"	0.0079 +/- 0.0020"	3.7 - 4.4 Fr
NFP017	2.00 mm	1.40 mm	0.20 +/- 0.05 mm	0.0787"	0.0551"	0.0079 +/- 0.0020"	4.6 - 5.7 Fr
NFP021	2.30 mm	1.70 mm	0.25 +/- 0.05 mm	0.0906"	0.0669"	0.0098 +/- 0.0020"	5.5 - 6.6 Fr
NFP026	2.84 mm	2.10 mm	0.25 +/- 0.05 mm	0.1118"	0.0827"	0.0098 +/- 0.0020"	6.7 - 8.3 Fr
NFP030	3.38 mm	2.50 mm	0.25 +/- 0.05 mm	0.1331"	0.0984"	0.0098 +/- 0.0020"	7.9 - 9.9 Fr

*Junkosha is manufactured and inspected in mm and values in inch are for reference only. Measurements of each dimension are performed using metric based measuring instruments.
*The inch values are calculated from the actual measured values (metric values) using the following formula, with one digit below the displayed digit rounded to the nearest value (numbers 5 and above are rounded up, and numbers below 5 are rounded down). (Inch value) = (Metric value) ÷ 25.4
*Peelable Heat Shrink Tubing meets the requirement of USP Plastic Class VI. Certificates available on request.